



**The Dolomite Centre Ltd**  
Unit 1, Anglian Business Park, Orchard Road,  
Royston, Hertfordshire, SG8 5TW, UK

**T:** +44 (0)1763 242491  
**F:** +44 (0)1763 246125  
**E:** sales@dolomite-microfluidics.com  
**W:** www.dolomite-microfluidics.com

**Dolomite Microfluidics – North America Office**  
Blacktrace Inc, 29 Albion Place  
Charlestown, MA 02129, USA

**C:** 617 803 6655  
**T:** 617 848 1121  
**F:** 617 500 0136  
**E:** salesus@dolomite-microfluidics.com

## INFORMATION SHEET

Part name	Miniature peristaltic pump (pack of 5) Miniature peristaltic pump (pack of 10)	Part number	3200054 3200055
-----------	---	-------------	--------------------

### Description

The miniature peristaltic pump is a compact pump that produces flow at 0.45 ml/min. It has very low power consumption of 0.12W and is supplied with silicone tubing fitted.



Miniature peristaltic pump

### Benefits

- Compact size
- Low power consumption
- Lightweight

	Specification	Value
1	Dimensions	30 x 12 x 14 mm
2	Flow Rate	0.45 ml/min
3	Discharge Pressure	0 – 50 kPa
4	Tube Size	OD 2.5 mm, ID 1.5 mm
5	Tube Material	Silicone
6	Motor	DC Geared Motor
7	Voltage	3V DC
8	Power Consumption	0.12 W
9	Weight	11 g

### Geometry

